

# Rapid methods for the analysis of environmental, decommissioning and bioassay samples – an update

S. Happel

# Actinides and Sr in soil, food, concrete and brick samples (decommissioning and emergency samples)

SL Maxwell, BK Culligan, A Kelsey-Wall, PJ Shaw: Rapid radiochemical method for determination of actinides in emergency concrete and brick samples. *Anal Chim Acta.*, 701(1):2011;112-118.

SL Maxwell, BK Culligan, A Kelsey-Wall, PJ Shaw: Rapid determination of actinides in emergency food samples, *J. Radioanal. Nucl. Chem.*, 292(1), 2011, 339-347

S. L. Maxwell and B.K. Culligan: Rapid Method for Determination of Actinides in Fecal Samples, 31/10/12, 58th Annual RRMCM, Fort Collins, CO October 29 to November 2, 2012

## **Rapid Method for Sodium Hydroxide Fusion of Concrete and Brick Matrices Prior to Americium, Plutonium, Strontium, Radium, and Uranium Analyses for Environmental Remediation Following Radiological Incidents**


**U.S. Environmental Protection Agency**

**Office of Air and Radiation  
Office of Radiation and Indoor Air  
National Analytical Radiation Environmental Laboratory  
Montgomery, AL 36115**

**Office of Research and Development  
National Homeland Security Research Center  
Cincinnati, OH 45268**



# Sample preparation




# Fluoride co-precipitation


- Carrier: Lanthanides and/or Ca
- Co-precipitate actinides and Sr
- U(IV) precipitates, U(VI) does not
  - Control of oxidation state allows for U discrimination
  - Ti(III) as reducing agent,  $H_2O_2$  as oxidizing agent
- Many matrix elements do not co-precipitate
  - Fe, Al, Ti,...
- Problematic: Use of HF
  - can be replaced by  $NaF/NH_4F$

# Separation scheme (Sr optional)

Pu, Np, Th: TEVA  
 Pu, Np, U, Th :TEVA+TRU  
 Pu, Np, Am, Cm :TEVA+DGA




# Rapid Purification of Pu for ICP-MS



**>10<sup>6</sup> removal  
of U-238 from Pu-239**

Add UTEVA + DGA Cartridge

Split Cartridges  
Discard UTEVA



Inductively Coupled Plasma (ICP) Mass Spectrometer

Health Physics Journal, August 2011  
– Vol 101 - Issue 2 - pp 180-186



# Method performance (MAPEP 18 samples)

- Good agreement (bias  $15\% \leq B \leq -15\%$ )
- High yields for actinides, good yields for Sr

| Sample Code              | Am yield (%) | Pu yield (%) | U yield (%) | Sr yield (%) |
|--------------------------|--------------|--------------|-------------|--------------|
| MAPEP-18 soil            | 96.2±6.33    | 102.2±10.5   | 84.0±5.64   | 60.0±2.8     |
| MAPEP-20                 | na           | na           | na          | 66.0 +/- 6.0 |
| 10g baby food            | 84.6±7.5     | 93.5±8.1     | 77.9±13.1   | na           |
| 10g apple                | 93.4±9.1     | 97.5±12.1    | 88.9±10.9   | na           |
| 10g squash               | 88.5±3.5     | 97.5±5.9     | 77.9±13.1   | na           |
| MAPEP-18 spiked concrete | 85.3±6.5     | 89.6±7.9     | 76.9±4.4    | na           |
| MAPEP-18 spiked brick    | 93.7±2.9     | 94.7±9.0     | 88.1±5.4    | na           |
| NRIP fecal               | 82.7±3.9     | 96.4±8.2     | 62.5±7.2    | na           |

S. Maxwell, 2010/11

- Results in < 1d – 2d
- Method can be adapted to ICP-MS



# Actinides in seawater

J Radioanal Nucl Chem (2014) 300:1175–1189

DOI 10.1007/s10967-014-3079-0

---

## **Rapid determination of actinides in seawater samples**

**Sherrod L. Maxwell · Brian K. Culligan ·  
Jay B. Hutchison · Robin C. Utsey ·  
Daniel R. McAlister**

# Sample preparation

- Applicable to 40 – 80L samples
  - can be simplified for 8L samples
  - Less La -> less DGA necessary
- Sample preparation in 4 – 8h
- Tracer yields between 85% and 95%

Seawater Sample  
Acidify to ~pH 2 with 12M HCl 1 mL per L)  
Add Tracers ( $^{242}\text{Pu}$ ,  $^{243}\text{Am}$ )

Add 50 mg Fe and 1 ml 10%  $\text{TiCl}_3$  per liter seawater

Add 12-20 mg La, Fe and  $\text{TiCl}_3$ , Mix, wait ~5 min. then add 2.5 mL 14.5  $\text{NH}_4\text{OH}$  per liter sample, mix, allow to settle, pour off supernatant to ~4- 6L (pH= 8.8-9.0)

Transfer to four 500 mL centrifuge tubes  
Centrifuge 10 minutes, discard supernatant.  
Rinse each ppt. with 100 ml water (pH 8.8-8.9).  
Centrifuge and discard supernatant.

80L: prepare two 40L aliquots with only 10 mg La in each. Combine after initial HCL/HF steps

Add 100 mL 1.5M HCl to one tube, redissolve and transfer to 2<sup>nd</sup> tube, rinse first tube 2-3 times with 20 mL 1.5M HCl, Cap and mix

80L: combine two 40L  $\text{LaF}_3$  ppts in 1.5M HCL in one 500 ml tube and repeat  $\text{LaF}_3$  ppt to combine

Add 15 mL 10%  $\text{TiCl}_3$ , mix  
Add 50 mL 28M HF, mix, and allow to sit ~10-15 min.  
Centrifuge and discard supernatant

Redissolve in 10 mL 3M  $\text{HNO}_3$ -0.25M Boric acid,  
9 mL 7M  $\text{HNO}_3$ , 12 mL 2M  $\text{Al}(\text{NO}_3)_3$


Valence Adjust:  
0.2 mL 1.5M sulfamic acid  
0.5 mg Fe (as iron nitrate)  
1.25 ml 1.5M ascorbic acid, wait 3 min.  
1 mL 3.5M  $\text{NaNO}_2$  then 2mL 15.8M  $\text{HNO}_3$

Column Load Solution

Centrifuge and wet-ash any organic material solids with 15.8M  $\text{HNO}_3$ , 3M  $\text{HNO}_3$ -0.25M Boric acid and 5mL 30 wt%  $\text{H}_2\text{O}_2$  on hotplate.

Redissolve in 5mL warm 3M  $\text{HNO}_3$ -0.25M Boric acid and add to column load solution.

# U, Pu/Np in 8L samples




\*ICP-MS: Elute Pu/Np with 15 ml 0.05M HCl-0.005M HF-0.02M Hydroxylamine Hydrochloride (HH), wet-ash  
or

Move Pu to DGA Resin via 1 ml UTEVA Resin with 15-20 ml 3M HNO<sub>3</sub>-0.1M AA-0.02M Fe<sup>2+</sup>, rinse with 5ml 8M HNO<sub>3</sub>, 20 ml 0.1M HNO<sub>3</sub> and elute Pu with 5-10ml 0.02M HCl-0.005M HF-0.02M HH

\*ICP-MS: elute U with 0.01M ammonium bioxalate or wet-ash with HNO<sub>3</sub>

# Pu/Np, Am/Cm in 40 - 80L samples



When 40L seawater aliquots are combined to get 80L, 4mL DGA Resin is needed to handle the amount of La added. DGA rinse volumes and final eluent are increased proportionally.

- Methods tested up to 80L
- Resin combination in function of radionuclides to be determined
- Tests on spiked rerel samples (various volumes and activity levels)

➤ TEVA/TRU method (1 – 8L)

- Pu-239:  $R_C = 89.3 \pm 9.7\%$ , Bias:  $-3.5 \pm 3.7\%$
- Np-237:  $R_C = 89.8 \pm 9.4\%$ , Bias:  $1.8 \pm 11.2\%$
- U-238:  $R_C = 94.7 \pm 6.3\%$ , Bias:  $-4.4 \pm 2.6\%$


➤ TEVA/DGA method (16 – 80L)

- Pu-239:  $R_C = 86.4 \pm 4.6\%$ , Bias:  $-0.3 \pm 8.0\%$
- Am-241:  $R_C = 94.0 \pm 3.6\%$ , Bias:  $-0.0 \pm 1.0\%$
- Cm-244:  $R_C = 94.0 \pm 6.3\%$ , Bias:  $-$

➤ TEVA method (10 – 20L)

- Pu-239:  $R_C = 91.1 \pm 5.9\%$ , Bias:  $-3$
- Np-237:  $R_C = 91.1 \pm 5.9\%$ , Bias:  $-1$

- Detection limit in function of sample volume for 16h count



# Sr-90 in seawater samples

J Radioanal Nucl Chem  
DOI 10.1007/s10967-014-3391-8

---

## Rapid determination of $^{90}\text{Sr}$ in seawater samples

Sherrod L. Maxwell • Brian K. Culligan •  
Jay B. Hutchison • Robin C. Utsey •  
Daniel R. McAlister

# Sample preparation

- Challenge: seawater contains high amounts of Sr and Ca
  - Sample volume for use of Sr resin limited
- Applicable to 40L samples
  - Sr-90 via Y-90
  - Separation on DGA
  - No Sr-89!
- Sample preparation in <8h

Seawater Sample  
Acidify to ~pH 2 with 12M HCl (1 mL per L)  
Add stable Y carrier (1mg)

Add 10 mg La to each aliquot. Add Fe (50 mg per L). Mix, add 14.5 M NH<sub>4</sub>OH (1.75mL per L), mix, allow to settle, pour off supernate to ~2- 4L (pH= 8.8-9.0)





Transfer to four 500 mL centrifuge tubes  
Centrifuge 10 minutes, discard supernate.  
Rinse each ppt. with 100 mL water (pH 8.8-8.9).  
Mix, centrifuge and discard supernate.

Add 100 mL 1.5M HCl to one tube, redissolve and transfer to 2<sup>nd</sup> tube, rinse first tube 2-3 times with 20 mL 1.5M HCl, Cap and mix

Add 75 mg Ca and 50 mL 28M HF, mix, and allow to sit ~15 min.  
Centrifuge and discard supernate.

Redissolve in 10 mL 3M HNO<sub>3</sub>-0.25M Boric acid,  
10 mL 15.8M HNO<sub>3</sub>, 10 mL 2M Al(NO<sub>3</sub>)<sub>3</sub>

Column Load Solution






- Sample preparation:  $\text{Fe}(\text{OH})_3$  and  $\text{CaF}_2$  co-precipitation
- Y separation on DGA, yield via ICP-MS
- Method applicable up to 40L
  - 80L when two 40L aliquots are combined
- Tested on spiked real samples
  - Varying volumes (10 – 40L) and activity levels (14.1 – 296 mBq/L)
  - $\text{RC}(\text{Y}) = 84.2 \pm 8.1\%$ , Bias:  $-1.4 \pm 2.9\%$
- Detection limits depend on counting time and sample volume
  - For GPC e.g.:


**MDA ( $\mu\text{Bq/L}$ ) for 40L Aliquot  
vs. Sample Count Time (min)**

Assumes a Detector Background of 1.5 cpm



**MDA ( $\mu\text{Bq/L}$ ) for 240 min Count Time  
vs. Sample Volume (L)**

Assumes a Detector Background of 1.5 cpm



# Sr-89/90 via SR/DGA

- Publications of Vajda et al and Maxwell et al.
- Stacked SR/DGA cartridges
- Load from 8M HNO<sub>3</sub>
- Sr retained and purified on SR resin
- Y retained and purified on DGA resin
- Both fraction analysed by Cerenkov counting
  - Sr fraction: Sr-89
  - Y fraction: Y-90 -> Sr-90
- Advantageous in case of high Sr-89/Sr-90 activity ratios
  - > Y-91



**Savannah River  
National Laboratory™**

OPERATED BY SAVANNAH RIVER NUCLEAR SOLUTIONS

We put science to work.™


# Rapid Method for $^{226}\text{Ra}$ in Urine Samples

**Sherrod L. Maxwell**  
Senior Fellow Scientist

*Radiobioassay and Radiochemical Measurements Conference 2013 – Rohnert Park, Ca*  
10/25/13

# Rapid Sample Preparation for Ra-226 in urine

*Less Ca, PO<sub>4</sub> instead of carbonate*



# Rapid Column Separation for Ra-226 in urine



- Improvement of existing method (Maxwell et al. 2011) for Ra-226 in urine
  - Use calcium phosphate to lower blank levels and reduce Ca (less cation resin)
  - Use DGA Resin to make method more rugged regarding evaporation steps
  - Combine cation resin elution with DGA Resin final purification of Ra-226
    - Stacked elution means only one evaporation step
    - Effective removal of U, Th, Po isotopes
  - High chemical yields
    - Spiked 100 mL urine samples (N = 6), two activity levels
    - $c_A(\text{Ra-226}) = 75.5 \text{ mBq/L}^{-1} \pm 6.1\%$  , bias = 3.9%,  $R_C = 92.8\% \pm 3.2\%$
    - $c_A(\text{Ra-226}) = 17.9 \text{ mBq/L}^{-1} \pm 4.5\%$  , bias = -2.7%,  $R_C = 98.0\% \pm 2.6\%$
    - U + Po decontamination > 500
  - < 3 hours with simultaneous sample preparation -> Ra-224
- Ba-133: No waiting for in-growth (but Ra-225 can be used if preferred)
- Can be adapted to smaller or larger urine aliquots as needed
  - Smaller aliquot if less urine available (spot urine sample)
  - Large aliquot if lower MDA needed



**Savannah River  
National Laboratory™**

OPERATED BY SAVANNAH RIVER NUCLEAR SOLUTIONS


We put science to work.™

# Rapid Determination of $^{210}\text{Po}$ in Water Samples

**Sherrod L. Maxwell**  
Senior Fellow Scientist


*Radiobioassay and Radiochemical Measurements Conference 2013 – Rohnert Park, Ca*  
10/24/13

# Rapid Sample Preparation Method for Po-210 in Water






# Rapid Column Separation Method for Po-210




\*DGA Resin - removes Th, U, Ra, Pu, Am isotopes




# Po Isotopes - Bismuth Phosphate (125 $\mu$ g Bi) -1L Groundwater Sample

---



# Rapid Sequential Separation Method for Actinides and $^{210}\text{Po}$ in Water



- Rapid method:
  - Calcium phosphate precipitation *~15 min. to 1 hour (depending on sample volume)*
  - Separation on DGA Resin *~1 hour*
  - Microprecipitation using bismuth phosphate *~15 min*
  - Sample ready for counting in 1.5 – 2.5h
  
- Tested on 200 mL to 1000 mL samples
- Analysis of spiked groundwater samples (N = 6)
  - 200 mL:  $c_A(\text{Po-210}) = 308.1 \text{ mBq/L}^{-1} \pm 4.8\%$  , bias = -2,4%,  $R_C = 87.4\% \pm 5,8\%$
  - 1000 mL:  $c_A(\text{Po-210}) = 61.5 \text{ mBq/L}^{-1} \pm 5.1\%$  , bias = -2,8%,  $R_C = 85.0\% \pm 8.2\%$
  
- Sequential analysis of Po and actinides possible
- Analysis of spiked groundwater sample (N = 6)
  - $c_A(\text{Po-210}) = 1660.4 \text{ mBq/L}^{-1} \pm 2.9\%$  , bias = 4.9%,  $R_C = 81.5\% \pm 2.6\%$
  - $c_A(\text{Pu-238}) = 381.1 \text{ mBq/L}^{-1} \pm 4.0\%$  , bias = 3.0%,  $R_C = 93.4\% \pm 6.8\%$
  - $c_A(\text{Am-241}) = 380.8 \text{ mBq/L}^{-1} \pm 2.9\%$  , bias = 2.9%,  $R_C = 100.2\% \pm 6.9\%$
  - $c_A(\text{Cm-244}) = 327.9 \text{ mBq/L}^{-1} \pm 3.7\%$  , bias = 0.1%,  $R_C = 100.2\% \pm 6.9\%$
  - $c_A(\text{U-238}) = 628.8 \text{ mBq/L}^{-1} \pm 3.7\%$  , bias = -4.3%,  $R_C = 96.6\% \pm 2.5\%$

# Thank you for your attention!



# TRISKEM

SHARING INNOVATION



TRISKEM INTERNATIONAL

Parc de Lormandière Bât. C - Rue Maryse Bastié - Campus de Ker Lann - 35170 Bruz - FRANCE

Tel +33 (0)2.99.05.00.09 - Fax +33 (0)2.99.05.07.27 - [www.triskem-international.com](http://www.triskem-international.com) - email : [contact@triskem.fr](mailto:contact@triskem.fr)