

Rapid methods for the determination of actinides and Sr in environmental samples

Scope

- Actinides and Sr in aqueous samples
- Actinides and Sr in soil, food, concrete and brick samples
- Determination of radiostrontium in large soil and sea water samples

Rapid methods


- Emergency situation
 - Rapid results, high sample throughput
 - Evaluation of situation, decision on measures to be taken
 - Usually low sample mass and short counting times
 - Routine analysis
 - High sample throughput
 - Sample volume / mass and counting time depending on detection limits to be obtained
- SRNL (Sherrod Maxwell) very active in developing new methods)


Determination of Sr, Pu, Am and U in water and urine samples

S.L. Maxwell: Rapid Analysis of Emergency Urine and Water Samples, Journal of Radioanalytical and Nuclear Chemistry , 275(3), 2008, 497 - 502

Eichrom method ACW17 VBS: Am, Np, Pu, Th, Cm, U and Sr in Water (with Vacuum Box System), Version 1.0, Eichrom Technologies LLC, October 2006


- Stacked TEVA, TRU and Sr cartridges
- Separation in < 6h (vacuumbox / cartridges)
 - Flow rates: 1 mL.min⁻¹ (load and elution), 3 mL.min⁻¹ (rinse)
- Results can be obtained in < 8h
 - incl. measurement, « emergency level » ≈ several Bq.L⁻¹
- 1 L water (pH 2) or mineralised urine
- Addition of internal standards and Sr-carrier (or Sr-85)
- Ca-Phosphat co-precipitation
- Dissolution in 3M HNO₃ / 1M Al(NO₃)₃
- Redox (Pu(IV)): Fe(II) / NaNO₂
- Load through all 3 cartridges
- Rinse with 3M HNO₃





TEVA Resin


Trialkyl, methylammonium
nitrate (or chloride)



R = C₈H₁₇ and C₁₀H₂₁

- Extractant: Aliquat 336®
- TEVA: **TE**t**ra**V**al**ent **A**ctinides
- Retention of Pu(IV), Th(IV), Np(IV)
- Sr, Am(III) and U(VI) not retained
- Th elution with 9M HCl
- Np/Pu co-elution with 0.1M HCl / 0.05M HF / 0.03M TiCl₃ or rongalite

Acid dependency of k' for various ions at 23°C.




Horwitz, et al. (HP195)

TRU Resin




- **Extractant:** CMPO / TBP
- **TRansUranium** elements
- Retention of Am(III) and U(VI)
- Am elution with 4M HCl
- Rinse with 4M HCl / 0.2M HF for Th removal
- U elution with 0.1M ammonium oxalate

Acid dependency of k' for various ions at 23-25°C.
TRU Resin



Horwitz, et al. (HP193)

Sr Resin




Diluent: 1-octanol

Dietz et al. 2004

Figures 2 and 3


Acid dependency of k' for various ions at 23-25°C.
Sr Resin



Horwitz, et al., (HP292)

Figures 4 and 5

Acid dependency of k' for various ions at 23-25°C.
Sr Resin



Horwitz (HP199)

- Extractant: 1.0M 4,4'(5')-di-t-butylcyclohexano 18-crown-6 in 1-octanol.
- Retention of Sr (and Pb)
- Rinse with 10 mL 8M HNO_3 , 5 mL 3M HNO_3 / 0.05M oxalic acid (in case of presence of tetravalent actinides) and 5 mL 8M HNO_3
- Sr elution with 10 mL 0.05M HNO_3 (Pb remains retained)
- Counting via GPC (alternatively LSC)

NRIP-2008 Water Analysis Results

Nuclide	Avg. Difference Reported vs NIST	Avg. Difference Longer Recounts
Pu-238	13 %	6.3%
Pu-240	- 2.3%	-4.5%
Am-241	9.6%	1%
U-238	-0.5%	-5.4%
U-234	9.0%	-6.7%
Sr-90	-14 %	N/A

Actinides: 45 minute count time / Recounts: 2 hour count time

Separation time 3 – 4h

S. Maxwell at 2007

Actinides and Sr in soil, food, concrete and brick samples


SL Maxwell, BK Culligan, A Kelsey-Wall, PJ Shaw: Rapid radiochemical method for determination of actinides in emergency concrete and brick samples. *Anal Chim Acta.*, 701(1):2011;112-118.

SL Maxwell, BK Culligan, A Kelsey-Wall, PJ Shaw: Rapid determination of actinides in emergency food samples, *J. Radioanal. Nucl. Chem.*, 292(1), 2011, 339-347

- Methods can be adjusted for larger sample masses
- Addition of internal standards and Sr carrier (or Sr-85)
- Mineralisation in furnace at 700°C
- NaOH fusion
- Two co-precipitations for matrix removal
 - $\text{Fe}(\text{OH})_3$ / Ca-Phosphate
 - LaF_3 under reducing conditions ($\text{TiCl}_3 \rightarrow \text{U}(\text{IV})$)
- Dissolution in 3M HNO_3 / 1M $\text{Al}(\text{NO}_3)_3$ / 0.25M boric acid
- Redox (Pu(IV)): Fe(II) / NaNO_2

- Vacuumbox system
- Stacked TEVA, TRU and DGA Kartuschen -> actinide retention
- Rinse with 3M HNO₃
- Separation of the cartridges (TEVA and TRU/DGA)
 - Th, Pu (Np) purification via TEVA
 - Am/U purification via TRU/DGA
- Microprecipitation
- Eluates from sample load and first rinse (all cartridges) united and evaporated to dryness
- Sr purification on 3 mL Sr Resin column or cartridges (2 mL + 1 mL)


Sample preparation



S. Maxwell et al. 2011




< 1 day

Separation scheme (Sr optional)



DGA Resin

Extractant



- DGA, Normal (*N,N,N',N'*-tetra-*n*-octyldiglycolamide), TODGA = DN
- High uptake of Actinides from 3M HNO_3
- Am / U separation with 0.1M HNO_3
- Am retained from 3M HCl
- Fe removal with 1M HNO_3 and La removal with 0.05M HNO_3
- Am elution with 0.25M HCl

Method performance (MAPEP 18 samples)

- Good agreement (bias $15\% \leq B \leq -15\%$)
- High yields for actinides, good yields for Sr


Sample Code	Am yield (%)	Pu yield (%)	U yield (%)	Sr yield (%)
MAPEP-18 soil	96.2±6.33	102.2±10.5	84.0±5.64	60.0±2.8
MAPEP-20	na	na	na	66.0 +/- 6.0
10g baby food	84.6±7.5	93.5±8.1	77.9±13.1	na
10g apple	93.4±9.1	97.5±12.1	88.9±10.9	na
10g squash	88.5±3.5	97.5±5.9	77.9±13.1	na
MAPEP-18 concrete	85.3±6.5	89.6±7.9	76.9±4.4	na
MAPEP-18 brick	93.7±2.9	94.7±9.0	88.1±5.4	na

S. Maxwell, 2010/11


- MAPEP: 15 – 30 mBq.g⁻¹ level
- Results in < 1d
- Method can be adapted to ICP-MS

Alpha spectra


Pu and Np spectra for concrete sample



Am and Cm spectra for concrete sample



U isotope spectra for concrete sample



Alpha spectra, concrete sample
Microprecipitation

S. Maxwell, 2011



Sr in large soil and sea water samples

Maxwell S L, Culligan B K, Shaw S J: Rapid determination of radiostrontium in large soil samples, Journal of Radioanalytical and Nuclear Chemistry (2012) DOI10.1007/s10967-012-1863-2

Maxwell S L, Culligan B K: Rapid determination of Radiostrontium in large soil samples, 31/10/12, 58th Annual RRMCM, Fort Collins, CO October 29 to November 2, 2012

Maxwell S L, Culligan B K, Utsey R C: Rapid method for the determination of Radiostrontium in sea water samples, 31/10/12, 58th Annual RRMCM, Fort Collins, CO October 29 to November 2, 2012

Rapid Determination of Radiostrontium in Large Soil Samples


S. L. Maxwell and B. K. Culligan
Savannah River National Laboratory
October 31, 2012




Rapid Sr-89, Sr-90 Acid Leach Method for Larger Soil Aliquots

50 g samples

For smaller samples (25 g)
also use of HF



Rapid Sr-89, Sr-90 Column Separation Method for Soil



For large soil aliquots:

- * 4 ml Sr Resin
- * 10mL 3M HNO₃-0.05M Oxalic Acid
- * 18 mL 0.05M HNO₃

Column Load Solution

Sr Resin* (3mL) cartridge

Rinse column;
15mL 8M HNO₃;
*5mL 3M HNO₃-
0.05M Oxalic Acid;
15mL 8M HNO₃

Sample matrix and ⁹⁰Y removal;
start time for ⁹⁰Y ingrowth after
final 8M HNO₃ rinse

Elute Sr
* 15mL 0.05M HNO₃

Evaporate on planchet;
weigh for gravimetric
yield

Count ^{89,90}Sr by Gas
Flow Proportional
Counter

Recount after ⁹⁰Y ingrowth
10 days later to determine
⁸⁹Sr and ⁹⁰Sr



UGM Roma 21/02/2013

P2/

22

Maxwell et al. 2012

Results spiked soils samples

- Leached 50g soil samples:
 - 5.92 mBq.g⁻¹ level
 - Yield: 94.0% (+/- 2.6%, N=7), Bias: 0.43%,
 - MDC: 0.41 mBq.g⁻¹ for 90 min count
 - 11.84 mBq.g⁻¹ level
 - Yield: 89.6% (+/- 2.7%, N=7), Bias: -2.51%,
 - MDC: 0,17 mBq.g⁻¹ for 8h count
 - 59.2 mBq.g⁻¹ level
 - Yield: 89.3% (+/- 5,3%, N=7), bias: -2,36%
- 25g HF digest
 - 11.84 mBq.g⁻¹ level
 - Yield: 73,0% (+/- 5,1%, N=7), Bias: 6,14%,
 - MDC: 0,17 mBq.g⁻¹ for 8h count
- All results corrected for 1.35 mBq.g⁻¹ Sr-90 found in unspiked soil
- High Pb samples: hold before measurement or pass through 1 mL DGA cartridge for Bi removal




We Put Science To Work

Rapid Method for Determination of Radiostrontium in Seawater Samples

S. L. Maxwell, B.K. Culligan, and Robin. C. Utsey

Savannah River National Laboratory


October 31, 2012




Radiostrontium in sea water

- Sea water: 7 – 8 mg Sr / L, 400 mg Ca / L
- ICP-MS for yield
- Preferably samples > 1L for low detection limits
- Preconcentration by coprecipitation
- 2 options for separation:
 - Sr-89/90: combined Sr/DGA resins
 - Sr-90 only: DGA resin
- Measurement via GPC (LSC or Cerenkov also possible)



Rapid Sr-89/90 Sample Preparation Method for Seawater



*2L: higher PO₄




Rapid Column Separation Method



Load at 1 drop/sec on vacuum box


Shows stacked option

If present, ⁹¹Y is a problem!
 If agree, confirmatory but better to allow ingrowth, then separate Y-90



Maxwell et al. 2012





Y-90 ingrowth option

➤ SRNL preferred option for high Sr-89 and low Sr-90

➤ Y-91 no problem

Personal communication
Sherrod Maxwell, 2013



Results Sr-89/90 option

- 1L spiked sea water sample (7,66 mg Sr.L⁻¹), 4 mL Sr resin:
 - 2h counting time
 - Measurement via Sr-90: Yield: 88,8% (+/- 5,9%, N = 11), Bias: 1,2%
 - Measurement via Y-90: Yield: 95,0% (+/- 1,6%, N = 11), Bias: 3,1%
 - Good correspondance
- 2L sea water sample (7,70 mg Sr.L⁻¹), 6 mL Sr resin:
 - 2h counting time
 - Yield: 81,9% (+/- 5,0%, N = 4), Bias: 4,2%
- Measurement via GPC
- MDAs:
 - 1L sea water (2 x 2 mL cartridges)
 - 2L sea water (3 x 2 mL cartridges)
 - MDAs: 9.1 mBq.L⁻¹ (2h count), 4,4 mBq.L⁻¹ (8h count), 3.0 mB.L⁻¹ (1000 min count)
 - 6L sea water (three 2L aliquots combined after purification)
 - MDAs: 1.5 mBq.L⁻¹ (8h count), 1.0 mB.L⁻¹ (1000 min count)


Sr-89/90 option

- Similar methods suggested for environmental water samples
 - Groska J, Molnar Z, Bokori E, Vajda N: Simultaneous determination of ^{89}Sr and ^{90}Sr : comparison of methods and calculation techniques, Journal of Radioanalytical and Nuclear Chemistry, March 2012, Volume 291 (3),707-715
 - T. O'Brien et al.: The rapid determination of Strontium-89 and Strontium-90 in Environmental Samples. Presented at the MARC IX conferences, Kailua-Kona, USA, 29/03/12
- Measurement by Cerenkov counting possible
 - Sr-89 and Y-90 via Cerenkov
 - Very low interference of Sr-90 on Sr-89
 - Advantageous in case of high Sr-89/90 activity ratios

Rapid Column Separation for ^{90}Sr (^{90}Y) – DGA Only Option



Avoids stable Sr limitation



Sr-90 (Y-90) DGA resin only option

- **1 to 10 liter method (DGA Resin only)**
 - 2 liter aliquot requires one 2 ml DGA Resin cartridge
 - *MDA with GFPC and 120 minute count = 9.1 mBq/L*
 - *MDA with GFPC and 480 minute count = 4.4 mBq/L*
 - *MDA with GFPC and 1000 minute count = 3.0 mBq/L*
 - 10 liter aliquot (5 x 2 liter aliquots combined after purification)
 - *MDA with GFPC and 480 minute count = 0.88 mBq/L*
 - *MDA with GFPC and 1000 minute count = 0.61 mBq/L*
 - *< 1mBq/L ⁹⁰Sr with 10L seawater aliquot and < 6 hour sample preparation*



TRISKEM

SHARING INNOVATION

TRISKEM INTERNATIONAL

Parc de Lormandière Bât. C - Rue Maryse Bastié - Campus de Ker Lann - 35170 Bruz - FRANCE

Tel +33 (0)2.99.05.00.09 - Fax +33 (0)2.99.05.07.27 - www.triskem-international.com - email : contact@triskem.fr