

CHARACTERISATION OF A TBP RESIN AND ITS APPLICATION TO THE SEPARATION OF SN

Dirks C.^a, E. Kovács-Széles^b, A. Bombard^a, S. Happel^{a*}

^aTriskem International, Bruz, FRANCE

^bCentre for Energy Research, Hungarian Academy of Sciences, Budapest, HUNGARY,

* corresponding author: shappel@triskem.fr

Introduction

TBP is a widely employed reagent in liquid-liquid extraction, especially in the extraction of actinides, one of the most prominent examples of its use being the Purex process. A TBP based extraction chromatographic resin has been characterized with respect to its U capacity and the weight distribution ratios (D_W) of U, Th, Pu, Np and numerous other cations in different concentrations of HNO_3 and HCl. Based on obtained data two methods for the separation of Sn have been developed: one for decommissioning samples (including a matrix removal step and special focus on Sn/Te separation - important for MS determination of Sn-126) and the other for the production of radionuclides (Sn separation from up to 1g of Cd).

Characterisation




Figure 1: Weight distribution ratios D_W of actinides and selected cations on TBP resin in various HCl and HNO_3 concentrations

- > Batch experiments
 - > measurement via ICP-MS or LSC
 - > D_W values = 60-110 for tetra- and hexavalent actinides in HNO_3
 - > High D_W values for Pu at elevated HCl concentrations
 - > Overall moderate D_W compared to TEVA/UTEVA/TRU/DGA Resins
 - > Easy strip of actinides at lower acid concentrations
 - > Interesting selectivities for Sn, Ga, Sb and noble metals in HCl
 - > Maximum uptake of the TBP resin was determined to be in the order of 75 mg $\text{U} \cdot \text{g}^{-1}$ resin in 8M HNO_3

Two step procedure for decommissioning samples

Sn separation from up to 1g Cd




Figure 2: Elution study, matrix removal via anion exchange, 10g 1X8 Resin (Eichrom A8, 200 – 400 mesh)




Figure 3: Elution study, Sn separation, Sn purification via 2mL TBP Resin columns




Figure 4: Elution study, Sn separation from 1g Cd, 2 mL column

- Matrix removal via anion exchange
 - > Direct or after co-precipitation for gross-matrix removal
 - > Load and first rinse from 6M HCl
 - > Further rinse with 1M HCl
 - > Removal of great part of matrix elements
 - > Sn elution in 50 mL 1 or 2M HNO_3
 - > Sn yield > 90%
 - > Trace Fe, Ga and Te still present, In follows Sn
- Sn purification via 2 mL TBP column
 - > Evaporation of Sn fraction and redissolution in 2M HCl for direct load onto TBP
 - > Rinse with 2M HCl
 - > Further removes matrix and In
 - > Rinse with 1M HCl
 - > Removes trace Ga, Fe and In
 - > Sn quantitatively eluted in 10mL 0.1M HCl
 - > Overall high decon. factor for Te
- Sn separation of Sn from 1g of Cd
 - > Quantification by ICP-MS
 - > Separation following Fig. 3
 - > Multi-element solution plus 1g of Cd (as chloride)
 - > No Cd determination in Load and Rinse possible (too high content)
 - > No Sn breakthrough during load
 - > Sn quantitatively eluted in 5mL 0.1M HCl
 - > Cd content in Sn fraction in low μg

Conclusions

- > TBP resin characterized with respect to D_W values and maximum U uptake
- > Good selectivity for Pu(IV), Sn, Ga and Sb in HCl
- > High potential for Sn separation/purification
- > Method development via elution studies
- > Sn separation methods developed
 - > Separation of Sn from decommissioning samples via two step procedure (AIX / TBP)
 - > Separation of Sn from large amounts of Cd
 - > Potential application to noble metal and Sb separation/purification