



# TrisKem International

Development of new extraction chromatographic materials for use in radioanalytical chemistry and isotope production

Steffen Happel

02/09/2022




6th-INCC 2022

6th International Nuclear Chemistry Congress



- Overview new resins/applications/R&D cooperations
- Examples of new resins and methods: radioanalysis
  - TK200 for actinide preconcentration and separation
  - TK221 for actinide separation
  - TK400/ZR Resin for Fe separation
- Examples of new resins and methods: radioisotope production
  - Ga-68 from Zn targets (ZR/TK200 & TK400/TK200)
  - Cu-61/4 from solid Ni targets (TK201/TBP)
  - Tb-161 from Gd targets (TK221 & TK211/2)
- Other new resins and ongoing projects



# Recently developed new products/applications


| Product         | Application                                           |
|-----------------|-------------------------------------------------------|
| TK100/1 Resin   | Sr, Pb direct separation, Ra separation               |
| TK200 Resin     | Actinides, Ga-67/8, Zn                                |
| TK201 Resin     | Cu-61/4, Tc, Re, Pu                                   |
| TK202 Resin     | Tc, Re                                                |
| TK211/2/3 Resin | Lanthanide separation (e.g. nca Lu-177, Tb-161,...)   |
| TK221 Resin     | Actinides, lanthanides                                |
| TK400 Resin     | Ga-67/8, Pa-230/1, Fe, Nb, Mo, Po                     |
| ZR Resin        | Zr-89, Ga-67/8, Ti-44/5, Ge-68, Fe                    |
| TK-TcScint      | Tc-99 direct LSC measurement (with Uni Barcelona)     |
| DGA Sheets      | Functionalized TLC paper for RN QC (with CVUT Prague) |

**Upcoming products:** **TK102 Resin** (Sr, Pb, Ra/Ba separation), **TK225 Resin** (lanthanide removal from acidic effluents), **TK300 Resin** (Cs and Rb separation), **impregnated membrane filters** (TK100, TK201,...)

# Recently developed new products/applications

| Product                | Application                                                |
|------------------------|------------------------------------------------------------|
| TK100/1 Resin          | Sr, Pb direct separation, Ra separation                    |
| <b>TK200 Resin</b>     | <b>Actinides, Ga-67/8, Zn</b>                              |
| <b>TK201 Resin</b>     | <b>Cu-61/4, Tc, Re, Pu</b>                                 |
| TK202 Resin            | Tc, Re                                                     |
| <b>TK211/2/3 Resin</b> | <b>Lanthanide separation (e.g. nca Lu-177, Tb-161,...)</b> |
| <b>TK221 Resin</b>     | <b>Actinides, lanthanides</b>                              |
| <b>TK400 Resin</b>     | <b>Ga-67/8, Pa-230/1, Fe, Nb, Mo, Po</b>                   |
| <b>ZR Resin</b>        | <b>Zr-89, Ga-67/8, Ti-44/5, Ge-68, Fe</b>                  |
| TK-TcScint             | Tc-99 direct LSC measurement (with Uni Barcelona)          |
| DGA Sheets             | Functionalized TLC paper for RN QC (with CVUT Prague)      |

**Upcoming products:** **TK102 Resin** (Sr, Pb, Ra/Ba separation), **TK225 Resin** (lanthanide removal from acidic effluents), **TK300 Resin** (Cs and Rb separation), **impregnated membrane filters** (TK100, TK201,...)




- Based on TOPO extractant
- High retention of actinides
- Applications:
  - Use for very efficient U removal from Pu (Wang et al.)
  - Determination of Tc-99 in water samples (Ni Yuan et al.)
  - Ga-68 production (in combination with ZR Resin)
  - Actinide separation from water samples
- Extracts actinides even at pH 1 - 2 (nitric acid)
  - Preconcentration and purification of selected actinides on same column
  - 'In the field'?


# Actinides on TK200 – HNO<sub>3</sub>

(all data N. Vajda et al)


### Dw of U(VI) and Th on TK200, mL/g




### Dw of Pu on TK200, mL/g



### Dw of Np(IV) on TK200, mL/g




### Dw of Am on TK200, mL/g




# Actinides on TK200 – HCl

(all data N. Vajda et al)


### Dw of U(VI) and Th on TK200, mL/g




### Dw of Pu(IV) on TK200, mL/g



### Dw of Np(IV) on TK200, mL/g




### Dw of Am on TK200, mL/g




# Actinides on TK200 – Application


Elution of Th from 1g TK200 resin cartridge



Elution of Pu from 1g TK200 resin cartridge




Elution of U from 1g TK200 resin cartridge



- Data by N. Vajda (RadAnal)
- Method development on-going
- Direct load of U, Th and Pu from acidified water samples (here 900 mL on 2 mL TK200)
- Sequential separation on TK200
- Automatisations & 'in the field' preconcentration'?




- Optimisation of DGA Resin
  - Contains TO-DGA / phosphine-oxide
  - higher load and more radiolysis stable inert support
- Main application: Lu, Tb concentration from high acid and elution in small volume of dilute HCl
- Higher U retention than DGA
- Higher Am retention than TRU
  - Potential interest for Actinide separation, particularly Am
- Cooperation with Nora Vajda
- Method development for water samples
  - TEVA/TK221 method
  - Ideally later also soil and decommissioning samples



Papp, I., Vajda, N. & Happel, S. *J Radioanal Nucl Chem* (2022). <https://doi.org/10.1007/s10967-022-08389-9>

# TK221 Resin



Th and Pu removed with TEVA.

Am and U separation on TK221 Am elution before U

**Table 3** Recovery of actinide tracers from spiked water samples

|                   | Actinides determination |                    |
|-------------------|-------------------------|--------------------|
|                   | Without Np separation   | With Np separation |
|                   | Yield %                 | Yield %            |
| TAP water         |                         |                    |
| <sup>230</sup> Th | 90 ± 8                  | 86 ± 7             |
| <sup>239</sup> Pu | 108 ± 7                 | 95 ± 7             |
| <sup>237</sup> Np | —                       | 91 ± 9             |
| <sup>241</sup> Am | 103 ± 7                 | 97 ± 6             |
| <sup>233</sup> U  | 103 ± 7                 | 70 ± 7             |
| SEA water         |                         |                    |
| <sup>230</sup> Th | 71 ± 7                  | 61 ± 6             |
| <sup>239</sup> Pu | 91 ± 7                  | 87 ± 6             |
| <sup>237</sup> Np | —                       | 93 ± 8             |
| <sup>241</sup> Am | 89 ± 7                  | 92 ± 6             |
| <sup>233</sup> U  | 88 ± 7                  | 59 ± 6             |

| Analyte           | Target values               |                         | Measured values        |                      |                               | Relative bias % | MARB <sup>a</sup> % | Z-score <sup>b</sup> | Test evaluation |
|-------------------|-----------------------------|-------------------------|------------------------|----------------------|-------------------------------|-----------------|---------------------|----------------------|-----------------|
|                   | Mean activity concentration | Standard deviation (sd) | Activity concentration | Standard uncertainty | Relative standard uncertainty |                 |                     |                      |                 |
|                   | Bq/kg                       | Bq/kg                   | Bq/kg                  | Bq/kg                | %                             |                 |                     |                      |                 |
| <sup>239</sup> Pu | 5.93                        | 2.27                    | 5.09                   | 0.24                 | 4.7                           | 14              | 25                  | 0.37                 | Accepted        |
| <sup>241</sup> Am | 4.85                        | 0.57                    | 4.73                   | 0.15                 | 3.2                           | 2.5             | 30                  | 0.21                 | Accepted        |
| <sup>244</sup> Cm | 7.02                        | 2                       | 7.19                   | 0.34                 | 4.7                           | 2.4             | 25                  | 0.09                 | Accepted        |



<sup>a</sup>Maximum Acceptable Relative Bias

<sup>b</sup>Z = |X<sub>reported</sub> - X<sub>target</sub>| / s<sub>dtarget</sub>



- Method tested on spiked tap and sea water samples
- High yields (88+ for U and Am)
- Analysis of IAEA-TEL-2021-03 WWOPT successful
- Next: use for solid samples?

# TK400 Resin

- Long chained alcohol
  - First work by Knight et al. on Np/Pa separation
- Retention only at high HCl concentration (>6M HCl), elution in low HCl, water,...
- Main application: Pa separation
- Also retains Mo, Fe, Po, Ga, Nb, Nb,...
- Higher Fe capacity than e.g. TRU Resin (~15mg Fe/g TK400)





# Fe/Mo/Nb separation



- Recovery of Fe/Nb/Mo from high HCl on TK400
- Majority of other elements removed during load and rinses (9M and 5M HCl)
- Fe/Nb and Mo eluted in dilute HCl => separation on ZR Resin
- Can also be used to remove Nb from Zr (e.g. stacked TK400/UTEVA) or Pu-241

# Ga-68 separation from Zn targets

- Irradiation of Zn-68 targets in cyclotron
- Ga-68 separation on ZR Resin
  - No selectivity for Zn (target material)
  - Loading possible from:
    - dilute acid (**liquid targets => typically HNO<sub>3</sub>**)
    - >6M HCl (**solid targets**)
  - Rinse under loading condition
  - Elution with ~1 - 2M HCl
  - Too acidic for injection or labelling



- Conversion necessary
  - Evaporation & dissolution difficult to automatize
- Easier => use of another resin
- TK200 Resin load from 1.5M HCl
- Rinse with 1.5M HCl
- Elution in 2 – 3 BV water, dilute acid,..

# Cyclotron production of Ga-68

Rodnick et al. *EJNMMI Radiopharmacy and Chemistry* (2020) 5:25  
<https://doi.org/10.1186/s41181-020-00106-9>

EJNMMI Radiopharmacy  
and Chemistry

RESEARCH ARTICLE

Open Access



## Cyclotron-based production of $^{68}\text{Ga}$ , $[^{68}\text{Ga}]\text{GaCl}_3$ , and $[^{68}\text{Ga}]\text{Ga-PSMA-11}$ from a liquid target

Melissa E. Rodnick<sup>1</sup>, Carina Sollert<sup>2</sup>, Daniela Stark<sup>3</sup>, Mara Clark<sup>1</sup>, Andrew Katsifis<sup>3</sup>, Brian G. Hockley<sup>1</sup>, D. Christian Parr<sup>2</sup>, Jens Frigell<sup>2</sup>, Bradford D. Henderson<sup>1</sup>, Monica Abghari-Gerst<sup>1</sup>, Morand R. Piert<sup>1</sup>, Michael J. Fulham<sup>4</sup>, Stefan Eberl<sup>5</sup>, Katherine Gagnon<sup>2</sup> and Peter J. H. Scott<sup>1\*</sup>




Fig. 4 Two-column approach for  $^{68}\text{Ga}$  chemical separation

**Table 1** High level schemes of  $[^{68}\text{Ga}]\text{GaCl}_3$  purifications

|                                  | Scheme A*                  | Scheme B                                                        |
|----------------------------------|----------------------------|-----------------------------------------------------------------|
| ① ZR Load                        | < 0.1 M $\text{HNO}_3$     |                                                                 |
| ② ZR Wash                        | 15 mL 0.1 M $\text{HNO}_3$ |                                                                 |
| ③ ZR Elution / Trapping on TK200 | 5–6 mL ~ 1.75 M HCl        |                                                                 |
| ④ TK Wash                        | –                          | 3.5 mL 2.0M NaCl in 0.13M HCl                                   |
| ⑤ TK Elution                     | $\text{H}_2\text{O}$       | 1–2 mL $\text{H}_2\text{O}$ followed by dilute HCl to formulate |

\*Process as reported previously (Nair et al. 2017)

- J. Kumlin et al.
- ZR, LN & TK200 for solid targets

ORIGINAL RESEARCH

### Multi-Curie Production of Gallium-68 on a Biomedical Cyclotron and Automated Radiolabelling of PSMA-11 and DOTATATE

> Helge Thisgaard, Joel Kumlin, Niels Langkjær, Jansen Chua, Brian Hook, Mikael Jensen, Amir Kassaian, Stefan Zeisler, Sogol Borjian, Michael Cross, Paul Schaffer, Johan Hygum Dam

DOI: 10.21203/rs.3.rs-70698/v1 [Download PDF](#)

- High Ga-68 activities
- ARTMS/Odense: 10 Ci production
- W. Tieu et al. use of single TK400 cartridge for solid Zn targets
- Svedjehed et al. use of TK400/A8/TK200 for solid Zn targets


Demystifying solid targets: Simple and rapid distribution-scale production of  $[^{68}\text{Ga}]\text{GaCl}_3$  and  $[^{68}\text{Ga}]\text{Ga-PSMA-11}$

Johan Svedjehed, Martin Pärnaste, Katherine Gagnon [✉](#)

# Cu-64 separation on TK201

- Cu-64 separation from solid Ni-64 targets
  - Original method development:
    - Target dissolution in high HCl
    - Load and rinse at 6M HCl
      - Ni removal and recovery/recycling
    - Co elution with 4 – 5M HCl
    - Cu elution with 0.5M HCl
      - Zn remains retained (Ga and Fe partially co-elute)

=> requires further treatment



- Improvements:
    - Preferred alternative: Use of TBP (or TK400) upfront for Fe/Ga removal
- => allows for Cu elution in 0.05M HCl

Svedjehed et al. *EJNMMI Radiopharmacy and Chemistry* (2020) 5:21  
<https://doi.org/10.1186/s41181-020-00108-7>

(2020) 5:21


*EJNMMI Radiopharmacy and Chemistry*

RESEARCH ARTICLE

Open Access

Automated, cassette-based isolation and formulation of high-purity [<sup>61</sup>Cu]CuCl<sub>2</sub> from solid Ni targets


Johan Svedjehed<sup>1</sup>, Christopher J. Kutyrail<sup>2</sup>, Jonathan W. Engle<sup>2,3</sup> and Katherine Gagnon<sup>1\*</sup>



- Gagnon et al. use of NaCl/HCl for better pH control of eluate

# Tb-161 separation from Gd targets

- nca Lu-177 still more frequently used but Tb-161 getting strong interest
  - Part of the ‘Swiss knife of nuclear medicine’ => Tb isotopes
- Similar production for both



| Tb 149        |               | Tb 152                     |                  | Tb 155       | Tb 161                 |
|---------------|---------------|----------------------------|------------------|--------------|------------------------|
| 4.2m          | 4.1 h         | 4.2m                       | 17.5 h           | 5.32 d       | 6.90 d                 |
| $\epsilon$    | $\epsilon$    | $\gamma$ 283;              | $\epsilon$       | $\epsilon$   | $\beta^-$ 0.5; 0.6...  |
| $\beta^+$     | $\alpha$ 3.97 | 160...                     | $\beta^+$ 2.8... | $\gamma$ 87; | $\gamma$ 26; 49; 75... |
| $\alpha$ 3.99 | $\beta^+$ 1.8 | $\epsilon$ ; $\beta^+$ ... | $\gamma$ 344;    | 105;...      | $\epsilon$             |
| $\gamma$ 796; | $\gamma$ 352; | $\gamma$ 344;              | 586;             | 180, 262     |                        |
| 165...        | 165...        | 411...                     | 271...           |              |                        |


**Terbium: a new ‘Swiss army knife’ for nuclear medicine**  
Source: <https://cerncourier.com/a/terbium-a-new-swiss-army-knife-for-nuclear-medicine/>

- Irradiation of several hundreds of mg or more
- Upscale on-going (incl. recycling) => typically 1g
- Use of prepacked PP columns
  - 4cm x 30cm (375 mL), 2.5cm x 30cm, 1.5cm x 30cm & 1.1cm x 30cm
  - Connection: ¼" 28G, up to ~10bar
  - HPLC pump
  - QC/CoA per column (peak asymmetry) for TK211/2/3







# Lanthanide separation on TK211/2/3



HDEHP (LN)




HEH[EHP] (LN2)



H[TMPeP] (LN3)



Extractants employed in **TK211/2/3**

- Extractant mixtures
- Optimized for high radiation stability




# Tb separation from 500 mg Gd targets

- Irradiated target typically oxide => dissolved in  $>3\text{M HNO}_3$
- Conversion via TK221 Resin
- Sequential separation on TK212/TK211
- Final conversion to dilute HCl on TK221 + trace nitrate removal on AIX




# Tb separation from 500 mg Gd targets

- Initial separation on TK212 – 150 mL column (30cm x 2.5cm)
- Large amount of Gd present leads to significant tailing
- Tb separation from Gd and Dy – ideally using online detection
- Fine purification on TK211 (29 mL)



Tb separation from 1000 mg Gd on TK212 (150 mL column)

# Tb purification on TK211



- Direct load of Tb fraction from TK212 onto TK211 (29 mL – 30cm x 1.1cm)
- Gd breakthrough during load & rinse with 0.5M HNO<sub>3</sub> (alternatively HCl)
- Tb elution: 2 options => 0.75M or 3.5M HNO<sub>3</sub> (3.5M HNO<sub>3</sub> preferable)
- Conversion to dilute HCl via TK221, A8 for nitrate removal
- Further improvement via EtOH addition?

- TK-TcScint
  - Direct Tc-99 measurement via impregnated PSm (=> Uni Barcelona)
- TK202 Resin
  - PEG based resin for Tc-99 separation from alkaline samples (e.g. after alkaline fustion)
- TK300 Resin
  - Cs and Rb selective resin based on a Calixarene
- TK102 Resin
  - Sr, Pb s& Ba/Ra separation
  - Higher  $k'$  and capacity than SR Resin
- Impregnated membrane filters
  - TK100, TK201, ...
  - In field preconcentration
  - DGT (Diffusive Gradients in Thin Films) => 'bio-availability'
- Radium
  - TK101
  - New resins and macrocycles
- Rapid tests => Test sticks
- Separation of DTM
  - SE Resin
  - Zr-93, Fe, Mo, Nb,...
- New Sheets (DGA, CU,...)
- Microfluidics
- ...

Thank you for your attention!





## SUBSCRIBE TO OUR NEWSLETTER

To keep updated with our latest developments, news and agenda for a year, subscribe to the TrisKem Infos [here](#)




# Tc-99 - new TK202 Resin

- Based on Polyethylene Glycol (PEG) grafted on inert support
- Aqueous Biphasic System (ABS)
- Retention of chaotropic anions like  $\text{TcO}_4^-$  in presence of kosmotropic anions ( $\text{SO}_4^{2-}$ ,  $\text{CO}_3^{2-}$ ,  $\text{OH}^-$ ,  $\text{MoO}_4^{2-}$ , ...)
  - Originally: Separation of Tc-99m from high masses of Mo




Dw values for Tc, Re and Mo on TK202 Resin, at varying NaOH concentrations. Tc data taken from Cieszykowska et al.(2).




Dw values for Tc in 5M NaOH using 40 mg TK202 Resin, increasing amounts of Mo. Data taken from Cieszykowska et al.

# Tc-99 - new TK202 Resin

- Tc retention from high NaOH (preferably 5 - 7M NaOH)
  - e.g. after alkaline fusion of decommissioning samples
- Re may serve as internal standard
- For high Mo samples:
  - Tc rec. > 90% for 6 – 8g Mo per g TK202
- Elution with water
  - Elution in small volume
    - Will still be alkaline
  - Pass through CEX for neutralisation and
  - through aluminium oxide for trace Mo removal and recovery as 0.9% NaCl solution




Re/Tc separation from Mo on TK202 Resin



Re separation from selected elements on 2 mL TK202 Resin cartridge, load and rinse at 1 BV/min, elution at 0.25 BV/min.





- Based on Plastic Scintillating microspheres (PSm)
- Scintillating beads impregnated with selective extractants
- Developed by García, Tarancón & Bagán
- Now available at TKI
- « TK **EI**Scint » range of products
- First: « TK **Tc**Scint »
  - Aliquot based > selectivity similar to TEVA
- Environmental/decommissioning monitoring => Tc-99 by LSC
- Direct measurement of cartridges after loading in LSC counter
  - Radioprotection, safety, hands-on-time, waste...
- Chemical yield via Re/ICP-MS in effluents







# Upcoming - TK300 Resin


- Macrocycle based Resin
- Cs and/or Rb separation
- Selectivity for Cs and Rb over other éléments tested in HNO<sub>3</sub> and HCl




Data provided by B. Russel (NPL)



# Upcoming - TK300 Resin



Elution study, Cs and Rb separation from selected elements on TK300 resin, loading from dilute acid.




Elution study, Cs and Rb separation from selected elements on TK300 resin, loading from 1M HNO<sub>3</sub>

- Separation of Cs and Rb
- Retention over wide pH range (up to 1M HNO<sub>3</sub>)
- Elution in >3M HNO<sub>3</sub>
- Cs/Rb separation possible
- Alternative => push resin into LSC vial (=>TEVA)
  - Discs?
- Drawbacks:
  - Limited Cs capacity
  - Interference by K
  - Limits use for environmental samples

➤ Rather suitable for decommissioning samples


- Range of **extractive membrane filters (MF)**
  - Rapid separation (up to 50 mL / min)
  - Preferably for use with water samples (0.5 – 5L)
  - Sampling in the field
  - Passive Sampling (DGT)
  - Under development:
    - TK201 (Tc)
    - TK100 (Sr, Pb), TK101 (Pb, Ra)
    - CL Resin (iodine)
    - Calixarenes (Ra)
    - TK300 (Cs)
    - ...
- Range of 'Test sticks'
  - Pieces of impregnated membrane on support
  - E.g. DGA for Ca,... (JCU, USouthampton)
  - Scintillating supports?



# Other separations on TK200


- Pd separation from Rh targets
- Main challenge: target dissolution & oxidation states
- Request: Pd separation from high  $H_2SO_4$
- Removal of  $H_2SO_4$  necessary
  - Rinse with 2M HCl
- Elution in acetate possible
  - To be optimized
- Separation on TK200 possible

## • Pt separation from Ir



- Pt/Ir separation. Elution study, ICP-MS measurement

## • Pd separation from Rh



- Pd/Rh separation. Elution study, ICP-MS measurement
- Pt separation from Ir targets
- Challenge oxidation state control
- Separation possible on TK200
- Alternative: use of TBP => Obata et al.
- $[^{188}, ^{189}, ^{191}\text{Pt}]$  cisplatin
- TBP and AIX based method
  - 3x 2 mL TBP cartridges followed by QMA cartridge

➤ Also Zn from Cu targets